



CanarE System Quick Start Guide



Read this guide carefully. If you have additional questions after reading this guide, refer to the Instructions for Use or contact Canary Medical.

1. System Components and Setup

1.1. CanarE System Components (Straps not shown):



1.2. CanarE Charging:

1. Connect the CanarE Charger to power with the supplied power cord and plug adapter. A USB power supply may also be used without the plug adapter.


NOTE: The Clinic Cart contains a power source for the CanarE Charger.

2. Place the CanarE flat surface down (logo face up) into the CanarE Charger holder. The CanarE contains an LED that will indicate its charge status while charging.
3. When ready to use the CanarE, remove it from the Charger. The CanarE automatically powers on.



LED Status	Definition
Constant Green	Charging: CanarE battery status is between 85%-100% charged
Flashing Green	Charging: CanarE battery status is between 40-99% charged
Flashing Orange	Charging: CanarE battery status is between 20-39% charged
Constant Orange	Charging: CanarE battery status is between 0-19% charged

1.3. Clinic Cart Setup:




2. Attaching the CanarE

1. Palpate the knee to locate the patient's left or right tibial tuberosity according to the Test Plan created with the Clinic Base Station System (Left Knee or Right Knee selected from the Clinic App Patient Test Plan).
2. Place the CanarE directly onto the patient's knee with the superior surface of the CanarE located on the inferior most aspect of the tibial tuberosity.

NOTE: Do NOT place the CanarE over patient clothing. This may cause an unintended shift of the CanarE during testing, which will affect test result data.

3. Align the CanarE alignment mark with the axis of the tibia by palpating the patient's tibia inferior of the CanarE.



4. Use the straps provided or self-adhesive tape to securely attach the CanarE to the patient:



- i. Loop one end of the strap through the lateral slot of the CanarE and secure the Velcro to the strap.



- ii. Wrap the other end of the strap around the patient's leg then through the other lateral slot of the CanarE.



- iii. Pull the strap tight and use the Velcro to secure it.

NOTE: Make sure the CanarE and strap are not too tight on the patient by asking the patient.



- iv. Check placement and adjust if needed.

NOTE: When pulling the strap tight to secure the CanarE to the patient, the CanarE may slightly shift on the patient's leg. Make sure that placement of the CanarE on the patient's leg has not changed. If the CanarE has shifted, readjust its position to the correct location. If the position is not corrected, this may affect the CanarE data collection.



Manufactured for:

Canary Medical USA LLC
2710 Loker Ave. West
Suite 350
Carlsbad, CA 92010
www.canarymedical.com
Customer Service: 1-833-722-6279 or 1-833-7CANARY

This material is intended for health care professionals. Distribution to any other recipient is prohibited.

© 2020-2022 Canary Medical Inc

PN. K03-CNE-300006 Rev. A